

Gel-U-Seat Lite General Use 2" Gel/Foam Wheelchair Cushion

HCPCS: E2601



FIGURE A



FIGURE B

- Offering great versatility, the Gel-U-Seat Lite provides comfort to almost all wheelchair users
- Unique dual-chambered gel bladder keeps the pelvis in contact with pressure redistributing cooling gel
- The cover is fluid resistant low-shear nylon zippered and removable with a durable vinyl bottom
- Secures easily to wheelchair with 12" tie straps
- Meets flammability standard TB117

PRODUCTS

Description

| | |
|--------|------------------------------------|
| 8040-1 | 2" Gel Cushion 16x16x2, 1 ea,4/CS |
| 8040-2 | 2" Gel Cushion 18x16x2, 1 ea,4/CS |
| 8040-3 | 2" Gel Cushion 20x16x2, 1 ea,4/CS |
| 8040-4 | 2" Gel Cushion 20x18x2, 1 ea,4/CS |
| 8040-5 | 2" Gel Cushion 18x18x2", 1 ea,4/CS |
| 8040-7 | 2" Gel Cushion 22x18x2", 1 ea/4CS |
| 8040-9 | 2" Gel Cushion 24x18x2", 1 ea,4/CS |

SPECIFICATIONS

Description

| | |
|-----------------|---|
| Cover Material | Fluid-Resistant Stretch Top Waterproof Vinyl Bottom |
| Warranty | 1 Year Limited |
| Weight Capacity | Up to 20" 275 lbs, 22" and up 500 lbs |