

## Gel-U-Seat Skin Protection 3" Gel/Foam Wheelchair Cushion

HCPCS: E2603



FIGURE A




FIGURE B

- The Gel-U-Seat skin protection cushion has a high density foam shell surrounding a uniquely designed gel bladder
- Dual chambered gel bladder provides proper immersion and envelopment for pressure redistribution
- Cover is water resistant and low shear stretch to reduce heat and moisture build-up
- Convenient 12" tie straps secure the cushion to the wheelchair
- Meets flammability standard TB117

### PRODUCTS

#### Description

|        |   |
|--------|---|
| 8046-1 | Gel-U-Seat W/C Cush, 16x16x3, 1ea,4/CS  |
| 8046-2 | Gel-U-Seat W/C Cush, 18x16x3, 1ea, 4/CS |
| 8046-3 | Gel-U-Seat W/C Cush, 20x16x3, 1ea, 4/CS |
| 8046-4 | Gel-U-Seat W/C Cush, 20x18x3, 1ea, 4/CS |
| 8046-5 | Gel-U-Seat W/C Cush, 22x18x3, 1ea,4/CS  |
| 8046-7 | Gel-U-Seat W/C Cush, 24x18x3, 1ea, 4/CS |

### SPECIFICATIONS

#### Description

|                 |   |
|-----------------|---|
| Cover Material  | Fluid Resistant Stretch Top   Waterproof Vinyl Bottom |
| Warranty        | 18 Months Limited                                     |
| Weight Capacity | Up to 20" 275 lbs, 22" and up 500 lbs                 |