

MEDNEB+ Compressor Nebulizer

HCPCS: E0570+A7003-1,E05070,A7005



FIGURE A

- Includes Airial compressor, disposable and optionable reusable nebulizer kit, 7' kink resistant tubing, mouthpiece, and extra filters
- Compact, Desk-Top Unit
- Piston-Driven System
- Optional carry bag with MQ5700B model

| PRODUCTS | Description |
|----------|--|
| MQ5600 | Compressor Nebulizer - MEDNEB 1/ea, 8/cv |
| MQ5700 | MEDNEB+ w/Reusable Neb Cup 1/ea, 8/cv |
| MQ5700B | MEDNEB+, Reusable Neb Cup 1/ea, 8/cv |

| SPECIFICATIONS | Description |
|-------------------------------------|-----------------------------|
| Compressor Free Air Flow | 5-8 LPM |
| Max Pressure | 39 psig or greater /269 KPa |
| Operating Humidity | 15-90% non-condensing |
| Operating Noise Level | <60 dBA |
| Operating Pressure Range | 8-13 psig or greater |
| Operating Temperature Range | 50 to 104°F (10 to 40°C) |
| Operating Voltage | 120VAC, single phase, 60Hz |
| Power Consumption | 180 VA |
| Size | 6.7" x 5.4" x 3.5" |
| Storage/Transport Humidity | 10-90% |
| Storage/Transport Temperature Range | -4 to 140°F (-20 to 60°C) |
| Warranty | Limited 5 Year |
| Weight | 2.6 lbs |

| ACCESSORIES | Description |
|-------------|-----------------------------------|
| 18254 | Reusable Neb Kit 10/case, CS |
| MASK 001A | Adult Aerosol mask CS/50, CS |
| NEB KIT 500 | Neb Kit,Inc Mouth pc&Tub50/cs, CS |