

## Model 177 Irrigation Syringe

HCPCS:



- The Model 177 Irrigation Syringe is efficiently used in various E.N.T irrigation procedures
- Small, malleable tube permits deep, yet painless introduction into external ear canal
- Regulate flow by placing finger over small hole in head of instrument if air pressure is 20 psi or under, then stream will stop when finger is removed
- Used by E.N.T doctors, medical doctors, research laboratories, veterinarians
- Can be operated by the No 630 or a No 1 bulb

| PRODUCTS     | Description                               |
|--------------|---|
| 177          | Irrigation Syringe                        |
| ACCESSORIES  | Description                               |
| 177-NKHR-605 | STRT FLW TIP ASSY RHODIUM PLATED 5 PK, EA |
| 115-QB       | 1-OZ CLEAR GLASS, EA                      |
| 177-OT-605   | OUTLET TUBE ASSEMBLY 5-PK, PKG            |
| 177-QB       | 6-OZ CLEAR GLASS BOTTLE, EA               |
| 114-S-605    | WASHER HERCULES 10-481 5-PK, PKG          |
| 630-605      | CUTOFF WITH TUBING 5-PK, PKG              |