

## Molded Foam General Use Wheelchair Cushion

HCPCS: E2601-1,E2601



FIGURE A

- A unique construction of high-resilient molded foam provides optimal support, stabilization and pressure redistribution
- Designed for low risk patients seeking comfort and support
- Removable, zippered, low-shear stretch cover is vapor permeable and water-resistant
- Cover has a "Non-Slip" vinyl bottom for comfort and safety
- Convenient 12" tie straps secure the cushion to wheelchair
- Meets flammability standard TB117

### PRODUCTS

#### Description

14880	W/C Cushion Molded, 16x16x2, 1/ea, 6/cv
14881	W/C Cushion Molded, 20x16x2, 1/ea, 6/cv
14887	W/C Cushion Molded, 18x16x2, 1/cv, 6/ct
14908	W/C Cushion Molded, 18x18x2, 1/ea, 6/cv
14909	W/C Cushion Molded, 20x18, 1/ea, 6/cv

### SPECIFICATIONS

#### Description

Cover Material	Fluid-Resistant Stretch Nylon
Warranty	18 Months Limited
Weight Capacity	250 lbs
Weight Capacity 18   20 (W)	250 lbs.   300 lbs.