

Nitro HD Hemi Height Rollator

HCPCS:



- 450 lb. weight capacity and low seat height
- Attractive, Euro-style design
- Brake cable inside frame for added safety
- Handle height easily adjusts with unique push button
- Unique suspension system for support and comfort
- Back support height easily adjusts with tool-free thumb screw
- Caster fork design enhances turning radius
- Large 10" front casters allow optimal steering and rolling comfort
- Easily folds with one hand to ultra-compact size for storage
- Under seat locking mechanism holds the Nitro closed when frame is folded
- Cross-brace design allows for side-to-side folding and added stability

SPECIFICATIONS

Description

Handle Height	35"–40"
Seat Dimensions	23" (W) x 10" (D)
Seat To Floor Height	19"
Warranty	Limited Lifetime 5 Years
Weight	20 lbs.
Weight Capacity	450 lbs