

## Oslo 3 Drawer Bedside Cabinet

HCPCS:



- High Pressure Thermo Fused Laminate
- Top High Pressure Laminate
- Locking top drawer is a standard
- Full Extension, Ball Bearing Drawer slides with rubber stops provide smooth operation, quiet closure
- Contemporary drawer handles with brushed nickel finish
- Ventilated back allows air flow throughout cabinet interior
- Available in three finish options (OSLO Cherry, OSLO Oak and OSLO Figured Mahogany)
- Specifications: 21"w x 18.00"d x 29.5"h

### PRODUCTS

### Description

NOBS3D-C	Oslo Bedside Cabinet, 3Drawer, Cherry
NOBS3D-FM	Oslo Bedside Cabinet, 3Drawer, FIG Mahogan
NOBS3D-O	Oslo Bedside Cabinet, 3Drawer, Oak