

Oxymizer® Disposable Oxygen Conserver

HCPCS:



FIGURE A

- P-224 and O-224 can be used with up to 15 lpm of continuous flow
- Facilitates the delivery of continuous high-flow oxygen therapy in a homecare, hospital, hospice or long-term care setting
- Delivers up to a 4:1 savings ratio
- Provides a comfortable alternative to a mask, allowing patients to eat, drink and talk
- Compatible with compressed gas cylinders, concentrators and liquid oxygen
- Proven efficacy for use at rest, during sleep and during exercise
- Disposable, does not require batteries
- Pendant model features a reservoir that can be worn underneath clothing

PRODUCTS

Description

O-224	Oxymizer Mustach- 24/Case
P-224	Oxymizer Pendant- 24/Case