

Pacifica 2 Nebulizer

HCPCS: E0570, A7003, A7013



- Ultra compact design makes it extremely portable
- Weighs only 3.3 lbs
- Powerful piston pump is efficient, dependable and quiet
- Maximum pressure: 30 psi
- Liter flow: 5-6 liters/minute
- 5 Year Limited Warranty
- Comes complete with 1 NEB KIT 500 packed in box

PRODUCTS	Description
18062	Pacifica 2 Nebulizer,1/ea

SPECIFICATIONS	Description
Electrical Requirements	120 Vac 60HZ
Litre Flow	5-6 liters/Minute
Maximum Pressure	30 PSI
Method of Nebulation	Piston
Operating Pressure	15 PSI
Particle Size	0.5 µm to 5 µm
Size	6.9" (W) x 5.1" (H) x 3.7" (D)
Warranty	Limited 5 Year
Weight	3.3 lbs (1.5Kg)

ACCESSORIES	Description
NEB KIT 500	Neb Kit,Inc Mouth pc&Tub50/cs, CS
MASK 001A	Adult Aerosol mask CS/50, CS
18254	Reusable Neb Kit 10/case, CS
18031	Nebulizer Carry Bag, 1/ea, EA