

## Penguin Pediatric Compressor Nebulizer

HCPCS: E0570 A7005-1, E0570, A7005



FIGURE A



FIGURE B

- Child friendly design
- Comes with disposable neb kit, pediatric mask, tubing, and a carry bag
- Durable, piston driven system

### PRODUCTS

### Description

MQ6002	Penguin Compressor Nebulizer 1/ea, 6/cv
MQ6002R	Penguin Neb,w/ReuseNebCup 1/ea, 6/cv
MQ6006	Penguin CompressorNebw/outBag 1/ea, 6/cv

### SPECIFICATIONS

### Description

Compressor Pressure Range	25 to 42 Psi
Liter Flow Range	5 ~ 8 lpm
Nebulizer Medication Capacity	5mL
Noise Level	55 dBA
Operating Humidity Range	10% to 95% RH
Operating Pressure Range	8 to 13 Psi
Operating Temperature Range	10° to 40°C (50° to 104°F)
Particle Size	0.5 to 10 µm
Power Consumption	80 VA
Power Requirements	AC 120V, 60Hz / 230V, 50Hz
Storage Humidity Range	10% to 95% RH
Storage Temperature Range	-20° to 70°C (-4° to 158° F)
Warranty	Limited 5 Year
Weight	3.52 lbs

### ACCESSORIES

### Description

18254	Reusable Neb Kit 10/case, CS
MASK 001P	Pediatric Aerosol mask CS/50, CS
NEB KIT 500	Neb Kit,Inc Mouth pc&Tub50/cs, CS