

## Premium Guard Water Mattress

HCPCS: E0198




FIGURE A



FIGURE B

- Three water baffles enhance movement and support to redistribute pressure and conform to varying body types
- Filling lines and "foot placement" lines provide for proper patient positioning
- Double-sealed seams provide maximum strength and durability
- Vinyl construction is durable and easy-to-clean with mild detergent
- Fill lines are conveniently located
- Fill hose sold separately (Figure A)

### PRODUCTS

#### Description

14400	Water Mattress,1/bx
-------	---------------------

### SPECIFICATIONS

#### Description

Warranty	Limited 90 Days on seams and valve
Weight Capacity	250 lbs
Weight Capacity (13046   13046-HD)	250 lbs.

### ACCESSORIES

#### Description

14401	Fill Hose for 14400, 1/ea, EA
-------	-------------------------------