

## RhinoGuard Atomizer Tip Cover

HCPCS:



- The RhinoGuard Atomizer Tip Cover is designed to prevent cross contamination and its associated risks
- Inexpensive and convenient to use, it fits the tips of DeVilbiss Atomizer models 15, 163, 180, and 286 our most popular models
- Precision-molded of soft, flexible polymer, the protectors help prevent trauma and seal the nasal passage during use
- Flared shield reduces splash-back for improved clinician protection
- Helps keep the atomizers tip in the correct location during use helping meet the challenge of treating pediatric and other active patients
- Enables more efficient administration of fine sprays

PRODUCTS	Description
RG250	RHINO GUARD, 250 PACK