

Deluxe Digital Electronic Muscle Stimulator

HCPCS:



- Three Stimulation Modes: Constant, Synchronous, and Alternate
- Timer: Adjustable 15, 30, 45, 60 minutes, or continuous
- Digital wave technology provides a highly efficient delivery of the wave signal
- Wave Form: Asymmetrical rectangular biphasic pulse
- Pulse frequency adjustable from 2-150 Hz
- Pulse Amplitude: Adjustable 0-100 mA into 500 ohm load each channel
- Pulse width is adjustable from 50 uS ~ 300 uS, 10 uS/step
- Comes complete with carrying case, 2 lead wires, 9-volt battery, 1 package of 4 pre-gelled electrodes and operating manual

SPECIFICATIONS

Description

Warranty

3 Year Limited