

## DeVilbiss Traveler® Portable Compressor Nebulizer System

HCPCS: E0570,A7005,A7013



FIGURE A



FIGURE B



FIGURE C

- Delivers a consistent and concentrated treatment of small particles that penetrate deeply into the lungs, delivering medication where its needed most
- Easy, one-touch operation
- Fast treatment time (approximately 10 minutes) with DeVilbiss 900 Series Reusable Nebulizer encourages patient compliance
- Can be used with any manufacturers nebulizer treatment times will vary
- 45-minute battery (NiMH) life far outlasts treatment time
- Includes DeVilbiss Traveler® compressor, disposable and reusable nebulizers, five extra filters, extra filter door, tubing adapter, adult and pediatric masks, AC adapter, DC car adapter, carrying case and instruction guide

### PRODUCTS

#### Description

6910D-DR	PORT.COMP. TRAVELER W/O BATT
6910P-DR	Traveler w/ Btty

### SPECIFICATIONS

#### Description

Battery Life	45 Minutes
Battery Warranty	1 Year on Battery
Compressor Warranty	3 Years
Dimensions	3.34" (W) x 2.2" (H) x 3.8" (D) Without Battery   3.34" (W) x 2.2" (H) x 5.3" (D) With Battery
Electrical	100-240V, 50-60Hz, 12VDC, 1.5A
Free Air Flow	6 lpm or Greater
MMAD	5 Microns
Max Pressure	30 psig or Greater
Sound Level	54 dBA
Weight	.8 lbs without battery   1.5 lbs with battery

### ACCESSORIES

#### Description

3655D-621	VixOne Disp. Nebulizer, CS 50, CS
6910D-601	REPLACEMENT BATTERY PACK, EA
DL1050	PEDIATRIC CHARACTER MASK (50-PK), CS
6910D-602	AC ADAPTER, EA
6910D-603	DC ACCESSORY POWER CORD, EA
6910D-605	INLET FILTERS (5) & COVER (1), EA