

PulmoMate® Compressor Nebulizer System

HCPCS: E0570



- The sleek case design of the PulmoMate® Compressor Nebulizer System fits easily into luggage, backpacks and diaper bags for easy transport. The DeVilbiss PulmoMate features an updated compressor/motor for long-term durability and performance.
- Robust quality and durability
- Consistent performance over the life of the unit
- Metal frame motor construction
- Integrated handle for easy carrying
- Large built-in compartment for storing nebulizers and tubing
- Intuitive design is easy-to-use and maintain

| PRODUCTS | Description |
|----------|------------------------------|
| 4650D | PulmoMate w/ Disp. Nebulizer |

| SPECIFICATIONS | Description |
|----------------|---------------------------------|
| Dimensions | 7.0" (W) x 3.8" (H) x 13.0" (D) |
| Free Air Flow | 13 lpm Typical |
| Max Pressure | 35 psig or Greater |
| Warranty | Limited 5 year |
| Weight | 5.3 lbs |

| ACCESSORIES | Description |
|-------------|--------------------------------------|
| 3655D-621 | VixOne Disp. Nebulizer, CS 50, CS |
| DL1050 | PEDIATRIC CHARACTER MASK (50-PK), CS |
| 5650D-602 | FILTER PKG - MM, 5/PK, PKG |
| 4650D-605 | STORAGE COMPARTMENT DOOR ASSY, EA |
| 5650D-612 | FILTER KIT(CAP&FILTER) (6/PK), PKG |
| 4650D-609 | PUSH-IN FEET (4-PK), PKG |
| 5650D-601 | FILTER PKG - DFT (25-PK), PKG |
| 4650D-611 | TUBE CONNECTOR W/PAL NUT 6-PKG, EA |