

## PulmoNeb® LT Compressor Nebulizer System

HCPCS: E0570,A7004



FIGURE A

- This affordable yet advanced system was developed to exceed market standards and provide superior clinical results and performance
- Small footprint and lightweight design make it easy-to-store and transport
- Darker body color maintains its appearance after a prolonged period of use
- Ball bearing connecting rod design equips the unit for years of dependable performance
- Internal wire guides eliminate interference with moving parts increasing reliability
- Top-of unit connector facilitates easy attachment/removal of tubing
- ETL listed as compliant to UL 60601-1 Medical Electrical Equipment Standard

PRODUCTS	Description
3655LT	PulmoNeb LT w/ Disp. Nebulizer,1/ea,6/CV
3655LTR	PulmoNeb LT w/ Disp. & Reus. Nebulizer

  

SPECIFICATIONS	Description
Dimensions	5.75" (W) x 4.0" (H) x 6.0" (D)
Free Air Flow	12.7 lpm Typical
Max Pressure	30 psig or Greater
Power Consumption	80 Watts Maximum
Sound Level	56-60 dBA
Warranty	Limited 5 Year
Weight	2.75 lbs

  

ACCESSORIES	Description
3655D-621	VixOne Disp. Nebulizer, CS 50, CS
DL1050	PEDIATRIC CHARACTER MASK (50-PK), CS
3655LT-601	LT FILTER, 5 PK, PKG
3655LT-602	LT FILTER, 25 PK, PKG
3655LTD-604	Snap In Feet Pkg 4, EA