

Bariatric Steel Forearm Crutch

HCPCS:



- Vinyl-coated, ergonomically contoured arm cuffs molded for comfort and stability
- Leg and forearm sections adjust independently for optimal sizing
- Vinyl hand grips are comfortable and durable
- Extra-large tips for added stability

PRODUCTS

Description

| PRODUCTS | Description |
|----------|------------------------------------|
| 10403HD | ForearmCrutch,Bariatric,Adult,1/cs |

SPECIFICATIONS

Description

| SPECIFICATIONS | Description |
|-----------------|------------------|
| Warranty | Limited Lifetime |
| Weight | 5.7 lbs. |
| Weight Capacity | 500 lbs |