

## Thigh Prompts

HCPCS: E1399-1,E1399



- Height adjustable
- Can be positioned further from or closer to the user to limit adduction/abduction
- Moves with users gait pattern allowing weight shifting from one leg to the other
- Prevents leg scissoring

### PRODUCTS

### Description

TK 1090 L	Thigh Prompts for Moxie, Lg,1pr/bag
TK 1090 S	Thigh Prompts for Moxie, Sm,1pr/bag

### SPECIFICATIONS

### Description

Distance From Clamp to Strap (TK 1090 S   TK 1090 L)	5"-8"   5"-11"
Maximum Length Circumference (TK 1090 S   TK 1090 L)	15"   20"