

## Trunk Support

HCPCS: E1399-1,E1399



FIGURE A

- Secures the torso at desired angle
- Height, depth, width and angle adjustable
- Padded for cushioned comfort
- Stays open for easy transfers
- Can be mounted towards the inside or outside of handlebar
- Adjustable straps provide a range of support from snug to loose depending on users capabilities
- In its lowest position, can be used as pelvic positioner
- Adjustments can be made while user is in the gait trainer
- Small trunk support should not be used on TK 3000

### PRODUCTS

### Description

TK 1080 L	Trunk Support for Moxie, Lg,1/bx
TK 1080 M	Trunk Support for Moxie, Med,1/bx
TK 1080 S	Trunk Support for Moxie, Sm,1/bx

### SPECIFICATIONS

### Description

Product Chest Circumference(TK 1080 S   TK 1080 M   TK 1080 L)	14"-28"   22"-40"   28"-50"
Product Height (TK 1080 S   TK 1080 M   TK 1080 L)	4.5"   6.5"   8"

### PARTS

### Description

TK 1082	CompAssemW/Brcket,Straps,Clips,TK 1080 S, EA
TK 1083	CompAssemW/Brcket,Straps,Clips,TK 1080 M, EA
TK 1084	CompAssemW/Brcket,Straps,Clips,TK 1080 L, EA
TK 1081	Height Adjustment Vertical Bar, TK 1080, EA