

## Hinged Toilet Seat Riser

HCPCS: E0244-1



FIGURE A



FIGURE B



FIGURE C

- Raises the height of the toilet 3.5" for individuals that require assistance sitting on and standing from the toilet
- Hinge design allows riser to be moved out of the way when not in use (Figure A)
- Heavy-duty molded plastic construction provides strength and durability
- Lightweight, portable, and easily attaches to toilet bowl
- Easy-to-install on most standard toilets
- Two models to accommodate either standard or elongated (Figure B) toilet bowls

### PRODUCTS

### Description

RTL12607	Standard Hinged Toilet Seat, 1 RTL/EA
----------	---------------------------------------

### SPECIFICATIONS

### Description

Seat Dimensions	13.5" (W) x 17.5" (D) x 3.5" (H)   14" (W) x 19.25" (D) x 3.5" (H)
Warranty	Limited Lifetime
Weight	5.07 lbs   5.62 lbs
Weight Capacity	250 lbs