

## Hospital Bed Linen Kit

HCPCS:



FIGURE A

- The Fitted Sheet is perfect for use with most standard manual, semi-electric and full electric beds
- Contains 2 fitted sheets that are ideal for standard hospital and homecare beds Fitted sheet measures 36" x 80" x 5" - 6"
- Sheets are made of a high quality cotton and polyester blend
- All components are machine washable

**PRODUCTS**

**Description**

PRODUCTS	Description
15030HBL	Hospital Bed Linen Kit,1 EA,10/CS
15030HBL-3684	Fitted Sheets, 36 x 84 x 8, 2/bx