

## Home Bed Assist Handle

HCPCS: E0315




FIGURE A



FIGURE B

- Assistance for getting in and out of bed
- Can be used on either side of the bed (Figure A)
- Removable handle
- Easy, quick, no tool assembly Comes in two pieces
- Safety strap wraps around mattress or box spring to ensure a secure fit
- Special extra-long safety strap enables user to be in a single, twin, double, queen or king size bed
- Made of powder coated steel for long life and easy cleaning
- Handle with multiple crossbars creates easy grip for any height use
- Retail Packaged

### PRODUCTS

### Description

15064	Bed Assist Handle,1/ea, RTL
-------	-----------------------------

### SPECIFICATIONS

### Description

Rail Dimensions	11.5" (W) x 20.75" (H) x 21.5" (D)
Warranty	Limited Lifetime
Weight Capacity	250 lbs