

T Style Half Rails

HCPCS: E0305



- Attractive brown-vein finish (Figure A)
- Constructed of 1" steel
- For use with "Pan-Style" beds and Bariatric Beds (item #s 15300, 15302 and 15303)
- Ideal protection without confining the patient
- New and improved design prevents rail from shaking
- Provides maximum protection when bed is elevated
- For use with most manufacturers' beds

PRODUCTS

Description

PRODUCTS	Description
15021BV	Half Rail,T-Style,BrownVein,1pr/cs