

## Deluxe Steel Drop-Arm Commode

HCPCS: E0165



- Easy-to-assemble frame
- Easy-to-release drop-arm mechanism allows for safe lateral patient transfers to and from commode
- Durable plastic snap-on seat and lid
- Easy-to-clean gray powder-coated steel finish
- Plastic armrests and back provide extra comfort and support
- Comes with 7 qt commode bucket with carry handle and splash shield

PRODUCTS	Description
11125KD-1	Deluxe Steel Drop-Arm Commode, CS
11125KD-2	Deluxe Steel Drop-Arm Commode, CS

  

SPECIFICATIONS	Description
Carton Shipping Weight	20   39 lbs.
Height	20" - 24"
Seat Dimensions	13.75" (W) x 16.5" (D)
Warranty	Limited Lifetime
Weight	18 lbs.
Weight Capacity	300 lbs
Width	19.75" between arms   23.5" Outside

  

ACCESSORIES	Description
RTL12085	CommodeLiners,EA 1BX of 12Liners,6BX/CS, CS