

## Steel Drop-Arm Commode with Wheels and Padded Armrests

HCPCS: E0165 OR E0166



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Easy-to-release arm mechanism allows for safe lateral patient transfers to and from commode
- Arm provides additional support while transferring
- Manufactured with 1" steel
- White powder coated frame is attractive and easy to maintain
- Padded arm rests provide extra comfort and support
- Comes with 7.5 qt. commode bucket with carry handle and splash shield
- Rust resistant, swivel casters (2 with locks)
- Limited Lifetime Warranty
- Made using recycled materials

### PRODUCTS

### Description

11101W-2	Commode,Drop Arm w/ Wheels,2/cs
----------	---------------------------------

### SPECIFICATIONS

### Description

Carton Shipping Weight	47 lbs.
Height	17" - 21"
Seat Dimensions	14.25" (W) x 16.75" (D)
Warranty	Limited Lifetime
Weight	22 lbs.
Weight Capacity	300 lbs
Width	18" between arms   22" Outside