

Airgo Comfort-Plus Cane with MiniQuad Ultra-stable Tip

HCPCS:



FIGURE A



FIGURE B

- Our popular offset Comfort-Plus cane now featuring our innovative MiniQuad ultra-stable tip
- MiniQuad tip is ultra-stable and supportive, while being 54% lighter than a large base quad cane
- The unique design re-centers itself for maximum stability
- Allows your cane to stand on its own
- The Comfort-Plus handle with a cushioned grip helps absorb shocks
- It is ideal for those with arthritis
- The 12 push-button height settings ensure this cane can be adjusted to fit you properly
- Built for longevity from ultra-durable lightweight anodized aluminum
- Includes a cane strap with reflective patch for added night time safety
- Fits most people 5' 6"
- AADL# C086

PRODUCTS

Description

730-857	AIRGO MINIQUAD CANE BLACK, 1EA, 6CV,RET
730-858	AIRGO MINIQUAD CANE BURG, 1EA, 6CV, RET
730-859	AIRGO MINIQUAD CANE BLUE, 1EA, 6CV, RET

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Weight Capacity	300 lbs
Weight Capacity (13046 13046-HD)	300 lbs (136 Kg)