

Steel Forearm Crutches

HCPCS: E0110 OR E0111



- Vinyl-coated, ergonomically contoured arm cuffs molded for comfort and stability
- Leg and forearm sections adjust independently for optimal sizing
- Vinyl hand grips are comfortable and durable
- Extra-large tips for added stability
- Ortho K Handle is hard plastic

PRODUCTS	Description
10403	Forearm Crutch, Adult, 1pr/bx
10405	Forearm Crutch, Tall 1pr/bx
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight	2.5 - 2.6 lbs.
Weight Capacity	300 lbs
ACCESSORIES	Description
10439-8	CrutchTips,10400,10406,10405Gry, 8/csRTL, CS