

## Steel Forearm Crutches



HCPCS:



FIGURE A



FIGURE B



FIGURE C

- Vinyl-coated, ergonomically contoured arm cuffs molded for comfort and stability
- Leg and forearm sections adjust independently for optimal sizing
- Vinyl hand grips are comfortable and durable
- Extra-large tips for added stability
- Ortho K Handle is hard plastic

### PRODUCTS

#### Description

10403	Forearm Crutch, Adult, 1pr/bx
10405	Forearm Crutch, Tall 1pr/bx

### SPECIFICATIONS

#### Description

Warranty	Limited Lifetime
Weight	2.5 - 2.6 lbs.
Weight Capacity	300 lbs

### ACCESSORIES

#### Description

10439-8	CrutchTips,10400,10406,10405Gry, 8/csRTL, CS
---------	--