

Airgo Claw Standing Cane Tip

HCPCS:



FIGURE A



- Lets your cane stand on its own! 360° of stability makes it ideal for folding canes
- 57% lighter than the MiniQuad weighing only 0.2 lb (92 grams)
- Feature a metal insert for durability
- Tool-free installation compatible with all canes with 3/4" (1,9 cm) cane tubing
- Even fits with Ice-O-Grips for winter safety

PRODUCTS

Description

735-002	Airgo Claw Tip, Black, 1EA,6CV,72CT,RET
---------	---

SPECIFICATIONS

Description

Warranty	Limited 90 Days
----------	-----------------