

## Respiratory Rollator Accessory Pack



HCPCS:



FIGURE A



FIGURE B

- Removable oxygen carry bag
- Multipurpose mount for cell phones, smart phones, and pulse oximeters
- Tubing and wire system to keep oxygen tubing and pulse oximeter wire probes from getting snagged
- Can fit up to two oxygen cylinders (ML6, M4, M9, M6) or one portable oxygen concentrator

PRODUCTS	Description
RESPROLL	Resp. Rollator Accessory Pack, 1/EA
RTL10555RD-RESP	Iwalker w/Resp Kit, Red

  

SPECIFICATIONS	Description
Dimensions	16" (W) x 6" (L) x 8.5" (H)
Handle Height	33.75" - 38.75"
Seat Dimensions	18" (W) x 9.75" (D)
Seat Height	20"
Unit Dimensions	27" (L) x 25.625" (W)
Warranty	1 Year Limited   Limited Lifetime   5 Years
Weight	1.4 lbs.   17 lbs.
Weight Capacity	300 lbs
Weight Capacity (13046   13046-HD)	300 lbs.
specifications	I-Walker rollator   accessory pack carry bag