

## Duet Rollator/Transport Chair, 8" Casters

HCPCS: E0143+E0156



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F

- Combines the features of a rollator and transport chair in one unit An individual can ambulate independently or be pushed safely by a caregiver
- To use as a transport chair the Duet has locking flip-up armrests (Figure A - Figure C) and a flip-down footrest (Figure D) to aid the user A safety lock release is used to store the arms and a hook and loop fastener strap is used to store the footrest while using the Duet for ambulation
- The tool-free two-position contoured padded backrest is reversible depending on whether you are using the Duet as a rollator or in transport chair mode
- Comes with padded seat with carry pouch for added privacy and security of personal items Inside the pouch is a handle to aid with the folding
- The folding frame has a lock-out feature to ensure the unit stays open while in use and ensures the unit will stay folded during transport
- 8" casters are ideal for indoor and outdoor use
- Handles are adjustable in height, with self-threading knob (Figure E)

### PRODUCTS

#### Description

795B	Duet Rollator/TransportChair,Blue,1/cs
795BK	Duet Rollator/TransportChair,Black,1/cs
795BU	Duet Rollator/TransportChair,Burgdy,1/cs

### SPECIFICATIONS

#### Description

Back Height	35"
Carton Shipping Weight	25 lbs.
Handle Height	31.5" - 37"
Overall Dimensions	25.5" (W) x 37" (H) x 26" (L)
Seat Dimensions	13.5" (W) x 13" (D)
Seat To Floor Height	21"
Seat Width	13.5"
Warranty	Limited Lifetime
Weight	20 lbs.
Weight Capacity	300 lbs

### PARTS

#### Description

9501A79512	Handlebar, Right, 795, 1/ea, EA
9502M7951L	Elbow Rest, Left, 795, 2M Ser.#, 1/ea, EA
9502M7951R	Elbow Rest, Right, 795, 2M Ser.#, 1/ea, EA
750FORK-2M	Front Fork, 795, 2M Ser.#, 1/ea, EA

PARTS	Description
9501A7953BK	FootrestAssm.Black,795 1A/2Y Ser.#,1ea, EA
9501A7953R-2M	FootrestAssembly,Red,795 2M/20hSN#, 1/ea, EA
9501A79511	Handlebar, Left, 795, 1/ea, EA
9501A7953B-2M	FootrestAssembly,Blue,7952M/20HSN#, 1/ea, EA
9501A79510-RP	Roll Pin Fork Fastener,795,1A/2Y SN#,1ea, EA
9501A79510-2M	Fork Stem Bolt, 795, 2M /20H SN#, 1/ea, EA
9501A79506	RearAxleBolt795,1A/2Y SN#,1/ea, EA
9501A79506-2M	Rear Axle Bolt, 795, 2M /20H SN#, 1/ea, EA
9501A79504	Tote 795 Rollator, 1A/2Y SN#,1ea, EA
750BCL-2M	HandBrakeAssy,Lt,795,2M/20H Ser.#, 1ea, EA
750BRAKE-2M	Universal hand brake assy, 795, 1/ea, EA
9500075005	Folding Hinge Bolt,750/795,1/ea, EA
9501A75013	Seat Bracket, 795 Rollator, 1/pr, EA
9502M75013	SeatBrackets&Screws,7952M/20H SN#, 1/pr, PR
9501A75008	3 Hole Wheel,8"Rear,750N,795,1/ea, EA
STDS1011	Wheel Bearing, 8mm x 22mm 1/ea, EA