

Travel Rollator, Black



HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Constructed of durable steel, the Travel Rollator is jam-packed with features to keep it a step ahead of the competition, including a 350 lb. weight capacity, small footprint for optimal shelf & storage space, & retail packaging to help your customers self-educate on product features & benefits. Introduce your customers to the perfect solution for maintaining their independence & mobility.
- Stay Independent & Mobile Wherever You Go
- Strong, durable steel frame supports up to 350 lbs
- Ergonomic hand grip with easy to use loop-lock handbrakes provides comfort and control
- Push down on the handbrake to engage parking brake feature
- For a customized fit, easily adjust the handle height with self-threading knob
- Padded backrest provides additional support when in a seated position, and can be folded up or down, or removed entirely
- Padded seat allows you to take a break when needed and remain comfortable during periods of extended sitting
- Convenient storage pouch under seat for easy access to your personal items
- 6" wheels are ideal for maneuvering both indoors and outdoors
- A compact and foldable design allows you to easily transport your rollator anywhere
- Quick and tool-free assembly

PRODUCTS	Description
RTL800KD-BK	RTL KD Rollator, 6", 4Whl Steel,Blk,1/ea

SPECIFICATIONS	Description
Casters	6"
Color	Black
Dimensions	22"W x 24"D
Handle Height	30" - 36"
Inside Hand Grip Width	18"
Material	Steel
Number of Wheels	4
Product Weight	17 lbs
Seat Depth	12"
Seat To Floor Height	20"
Seat Width	14"
Shipping Weight	19 lbs
Warranty	Limited Lifetime
Weight Capacity	350 lbs

ACCESSORIES	Description
STDS1040S	Cup Holder, Swivel,1BX,6/CS, BX

ACCESSORIES**Description**

RTL10254BFL	Univ Mobility Bag, Med, Blue Flrl, 1EA, EA
RTL10254GR	Univ Mobility Bag, Med, Army Green, 1EA, EA
RTL10254TFL	Univ Mobility Bag, Med, Trop Flrl, 1EA, EA
RTL10255NV	Univ Mobility Bag, Small, Navy, 1EA, EA
RTL10255PK	Univ Mobility Bag, Small, Pink, 1EA, EA
RTL10256NV	Univ MobilityBag, Med Flat, Navy, 1EA, EA
RTL10256PK	Univ MobilityBag, Med Flat, Pink, 1EA, EA
RTL10256TL	Univ MobilityBag, Med Flat, Teal, 1EA, EA