

## Bariatric Aluminum Folding Walker- Two Button

HCPCS: E0148



- For Bariatric use
- Steel legs and side braces
- Precision design provides additional strength while adding minimal weight
- Wider and deeper frame design accommodates individuals up to 500 lbs
- Easy push-button mechanisms may be operated by fingers, palms or side of hand
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Vinyl-contoured hand grip

### PRODUCTS

### Description

10220-1	Bariatric Aluminum Folding Walker- Two Button, EA
---------	---

### SPECIFICATIONS

### Description

Carton Shipping Weight	12   21   18 lbs.
Overall Width	25.25"   24.5"
Warranty	Limited Lifetime
Weight	9   8 lbs.
Weight Capacity	500 lbs
Width	20.5" inside hand grip   23" inside back legs