

Competitive Edge 2-Button Walker with 5" Wheel

HCPCS:



FIGURE A

- Easy push-button release mechanisms
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Sturdy 1" diameter construction ensures maximum strength
- Newly designed rear glide cap allows for walker to slide easily and smoothly over most surfaces
- Vinyl-contoured hand grip
- Comes with wheels (Replacement wheels item # 10109)

SPECIFICATIONS	Description
Carton Shipping Dimensions	10.6" (W) x 20" (L) x 31.2" (H)
Handle Height	32" 39"
Inside Hand Grip Width	17.5"
Product Weight	8 lbs.
Warranty	Limited Lifetime
Weight Capacity	350 lbs
Weight Capacity (13046 13046-HD)	350 lbs.
Width inside back legs	19.25"