

## Deluxe Folding Walker, Two Button

HCPCS: E0135



FIGURE A

- Easy push-button mechanisms may be operated by fingers, palms or side of hand
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Unique U-shaped frame design creates greater clearance
- The patient is able to bring the walker closer to assist in standing because of the improved crossframe (Figure A)
- Sturdy 1" diameter, aluminum construction ensures maximum strength while remaining lightweight

### PRODUCTS

#### Description

10200-1	Walker,2Btn.Adult,1/ea
10200-4	Deluxe Two Button Walker, 4/CS
10201-4	Deluxe Two Button Junior, 4/cs

### SPECIFICATIONS

#### Description

Adjustable Height	32" - 39"
Carton Shipping Weight	26   24 lbs.
Overall Width	22"
Warranty	Limited Lifetime
Weight	6 lbs.
Weight Capacity	350 lbs
Width	17" inside hand grip   19" inside back legs