

## Folding Walker, Two Button with 5" Wheels

HCPCS: E0143



FIGURE A

- Easy push-button release mechanisms
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Sturdy 1" diameter construction ensures maximum strength
- Newly designed rear glide cap allows for walker to slide easily and smoothly over most surfaces (Figure A)
- Vinyl-contoured hand grip
- Comes with wheels (Replacement wheels item # 10109)
- Limited Lifetime Warranty

PRODUCTS	Description
10244-4	Walker,Steel, Adult 2Button 5inWhls,4/CS

  

SPECIFICATIONS	Description
Carton Shipping Weight	9 lbs   33 lbs
Depth (opened at base)	17.5"
Height	32"-39"
Overall Width	24"
Warranty	Limited Lifetime
Weight	8 lbs
Weight Capacity	350 lbs
Width	17" inside hand grip   20" inside back legs

  

PARTS	Description
10109	Walker Wheels,Univrsl 5",1pr/bx, RTL, BX