

Airgo Comfort-Plus XC Premium Transport Chair

HCPCS:



- Designed to be as pleasurable as it is practical
- Angled handlebars provide more room for ones back, while caregivers will appreciate the comfortable ergo-form handgrips with convenient lockable hand brakes for increased safety
- Unique footrests can be conveniently locked to the side of the frame for easy storage
- Detachable flip back armrests also make it easier to pull up to a table or desk
- Heel-loops and seat belt included as standard features for increased safety
- Large 12" rear wheels with premium bearings for smoother rolling and increased comfort
- Ergo-form handgrips comfortably fit your palm and provide natural brake access
- Handgrip mounted brake levers for increased control
- Anatomically correct seat-to-armrest height for a naturally erect posture
- Detachable full flip-back armrests for effortless transfers
- Fold down back for more compact transport and storage
- Convenient backrest pouch for discreet storage

PRODUCTS

Description

700-830	Airgo Comfort-Plus XC TC, Black, 1 EA
---------	---------------------------------------

SPECIFICATIONS

Description

Carton Shipping Weight	23 lbs
Caster Size	8" and 12"
Product Weight	21 lbs
Seat Depth	16"
Seat To Floor Height	19"
Seat Width	19"
Warranty	Limited lifetime on frame 90 days on wearable components
Weight Capacity	300 lbs
Width	Open: 24" (W) x 33" (D) Closed: 10.5" (W) x 33" (D)