

Bariatric Steel Transport Chair

HCPCS: E1039



FIGURE A



FIGURE B

- Supports individuals up to 450 lbs
- Reinforced steel frame provides added support
- Dual, reinforced steel cross braces
- Available in blue or red frames
- 12" rear flat free wheels
- Comes with swing-away footrests
- Padded armrests provide additional comfort
- Heavy-duty, nylon reinforced upholstery with a back carry pouch
- Removable, reversible desk-length arms
- Companion activated wheel locks

PRODUCTS

Description

BTR20-B	TransportChr.HD,20"Blue,1/cs
BTR20-R	TransportChr.HD,20"Red,1/cs
BTR22-B	TransportChr.HD,Blue,22",1/cs
BTR22-R	TransportChr.HD,Red,22",1/cs

SPECIFICATIONS

Description

Back Height	16"
Carton Shipping Weight	50 lbs 54 lbs
Overall Height	37"
Overall Length w/ Riggings	38"
Overall Width	24.5" (open) 9" (closed) 26.5" (open) 11.5" (closed)
Seat Dimensions	20" (W) x 18" (D) 22" (W) x 18" (D)
Seat To Floor Height	19.5"
Warranty	Limited Lifetime
Weight	45 lbs 49 lbs
Weight Capacity	450 lbs

PARTS

Description

STDS1072-8	Seat Upholstery Screws, 5mm, 8/pkg, EA
STDS2K3003	HAND GRIP,BTR, 1pr, PR
STDS2K303R	Brake Assembly, BTR Right, 1/ea, EA
STDS2K303L	Brake Assembly, BTR Left, 1/ea, EA
STDS2K3002	REAR HUB CAP,BTR,6/bag, BAG
STDS2K3000	Rear Wheel for BTR22-B,12" Wheel,1/ea, EA

