

## Multi-Use Accessory Bag

HCPCS:



- Back strap attachment with secure, adjustable clip
- Securely attaches to wide range of mobility aids
- Spacious interior pocket with zipper closure and key ring
- Front pocket with magnetic closure for added security of belongings
- Reinforced back foundation maintains bag shape for added support
- Interior small pockets and zipper pocket
- Fully-lined pockets and interior
- Removable, optional cross-body carrying strap
- Treated with water resistant technology
- Spot clean with wet cloth
- 

PRODUCTS	Description
RTL10255NV	Univ Mobility Bag, Small, Navy, 1EA
RTL10255PK	Univ Mobility Bag, Small, Pink, 1EA

SPECIFICATIONS	Description
Dimensions	6.5" H x 10" W x 3" D
Retail Box Dimensions	6.1" L x 1.4" W x 11" H
Weight	.7 lbs.
Weight Capacity	3 lbs.