

Low Bed Overbed Table

HCPCS: E0274



- Designed specifically for use with low beds with 3.25" frame to floor clearance
- Walnut woodgrain melamine top
- Top can be raised or lowered in infinite settings between 19.5"- 27.5"
- Tabletop can be raised with the slightest upward pressure
- The frame is constructed with tubular steel and durable, attractive chrome finish
- Swivel casters allow for easy maneuverability

PRODUCTS	Description
13081	Overbed Table,Low Style,1/ea
SPECIFICATIONS	Description
Base Dimensions	15" x 28"
Carton Shipping Weight	28 lbs.
Height	19.5" - 27.5"
Tabletop Dimensions	15" x 30"
Warranty	Limited Lifetime
Weight	25 lbs.
Weight Capacity	100 lbs