

Hydraulic Deluxe Chrome-Plated Patient Lift

HCPCS: E0630



- Chrome-plated steel provides maximum strength
- Safely raises or lowers individuals from any stationary position
- 6 point Swivel Bar that can accommodate, 2 or 4 sling straps, plus 2 chain connections
- High-performance hydraulics raise or lower individuals gradually and safely
- Easy-to-operate caster brakes provide security
- Includes sling chains

| PRODUCTS | Description |
|----------|-------------------------------|
| 13023 | Patient Lift, New Style, 1/cs |

| SPECIFICATIONS | Description |
|------------------------|---|
| Base Height | 6" |
| Base Width | 22.5" - 42" |
| Boom Height | 17.5" - 62" |
| Carton Shipping Weight | 81 lbs |
| Casters | 5" |
| Warranty | Limited Lifetime (frame), 1 Year Limited (pump) |
| Weight | 72 lbs |
| Weight Capacity | 450 lbs |