

Patient Lift Scale

HCPCS:



FIGURE A

- Recommended for use on DRIVE LEVANTAR (FLNP500/FLP500) or GRAVIS (FLNP600/FLP600) floor lifts
- Style and color compliment LEVANTAR and GRAVIS floor lifts
- 660 lbs/300 kg max weight capacity
- Easy to install
- Will not fit on old style 13240,13242,13244 and 13245 DRIVE Floor lifts

PRODUCTS

Description

FLSCALE	Patient Lift Scale 660lbs capacity
---------	------------------------------------

SPECIFICATIONS

Description

Battery Life	Approximately 150 hours
Battery Type	3-AA batteries
Display Dimensions	3.2 in x 1.8 in x 5.5 in (81.3mm x 45.7 mm x 139.7 mm)
Power Supply	3-AA batteries
Size	3.2 in x 1.8 in x 5.5 in
Weight	1.37 lbs / 620 g
Range	600lbs / 300kg