


Digital Scale

HCPCS:



- Determines the weight of a person suspended from a portable patient lift mechanism
- Easy to set up and operate
- Operates with 1.9V battery
- Battery power is saved by automatic shut-off feature
- Can be used with most manufacturers lifts
- 13046 works with Drive's lifts: 13010SV, 13023, 13240, 13242
- 13046-HD works with Drive's lifts: 13244, 13245

PRODUCTS	Description
13046-HD	Digital Scale,HD 600 LB Lift,1/ea

SPECIFICATIONS	Description
Warranty	90 days
Weight Capacity	400 lbs
Weight Capacity (13046 13046-HD)	400 lbs 600 lbs.