

Shower Stand

HCPCS: E1399-1,E1399



- Sturdy aluminum stand
- Swivel locking casters
- Quick-release locks allow chair to be easily attached or removed
- Assembly required

PRODUCTS	Description
OT 8020	Otter ShowerStand,wheeled,1/cs

SPECIFICATIONS	Description
Carton Shipping Weight	20 lbs
Unit Dimensions	27" (W) x 21" (H) x 33" (L)
Weight	17 lbs
Weight Capacity	200 lbs

PARTS	Description
OT 0023	Release Lever Assembly, OT 8020,1/ea, EA
OT 8024	Replacement Screws for OT 8020, 4/bag, BAG
OT 8022	Caster Frame Insert, OT 8020, 1/ea, EA
OT 0022	Release Lever, OT8010/OT8020, 1/ea, EA
OT TIPS	Replacement Tips,Otter,1 set/bx, BX
OT 8023	EndCapOtterShowerStand,1ea, EA
OT 8021	Locking Caster for OT 8020,1/ea, EA