

Thigh Prompts

HCPCS: E1399-1, E1399



- Height adjustable
- Can be positioned further from or closer to the user to limit adduction/abduction
- Moves with users gait pattern allowing weight shifting from one leg to the other
- Prevents leg scissoring

PRODUCTS

Description

TK 1090 L	Thigh Prompts for Moxie, Lg, 1pr/bag
TK 1090 S	Thigh Prompts for Moxie, Sm, 1pr/bag

SPECIFICATIONS

Description

Distance From Clamp to Strap (TK 1090 S TK 1090 L)	5"-8" 5"-11"
Maximum Length Circumference (TK 1090 S TK 1090 L)	15" 20"