

Trunk Support

HCPCS: E1399-1,E1399



FIGURE A

- Secures the torso at desired angle
- Height, depth, width and angle adjustable
- Padded for cushioned comfort
- Stays open for easy transfers
- Can be mounted towards the inside or outside of handlebar
- Adjustable straps provide a range of support from snug to loose depending on users capabilities
- In its lowest position, can be used as pelvic positioner
- Adjustments can be made while user is in the gait trainer
- Small trunk support should not be used on TK 3000

PRODUCTS

Description

| | |
|-----------|-----------------------------------|
| TK 1080 L | Trunk Support for Moxie, Lg,1/bx |
| TK 1080 M | Trunk Support for Moxie, Med,1/bx |
| TK 1080 S | Trunk Support for Moxie, Sm,1/bx |

SPECIFICATIONS

Description

| | |
|--|-----------------------------|
| Product Chest Circumference(TK 1080 S TK 1080 M TK 1080 L) | 14"-28" 22"-40" 28"-50" |
| Product Height (TK 1080 S TK 1080 M TK 1080 L) | 4.5" 6.5" 8" |

PARTS

Description

| | |
|---------|--|
| TK 1082 | CompAssemW/Brcket,Straps,Clips,TK 1080 S, EA |
| TK 1083 | CompAssemW/Brcket,Straps,Clips,TK 1080 M, EA |
| TK 1084 | CompAssemW/Brcket,Straps,Clips,TK 1080 L, EA |
| TK 1081 | Height Adjustment Vertical Bar, TK 1080, EA |