

## Cast Protectors

HCPCS:



FIGURE A



FIGURE B




FIGURE C

- Keeps casts and bandages dry while bathing or showering
- Stretch tight fit eliminates the need for pumps, straps and adhesive tapes
- Reusable design folds compactly
- Latex Free

### PRODUCTS

### Description

RTLPC23401	Cast Protector, 20" Arm, 1 PC/RTL, 50/CV
RTLPC23402	Cast Protector, 24.5" Leg, 1 PC/RTL, 50/CV

### SPECIFICATIONS

### Description

Case Weight	6.6 lbs   9.15 lbs
Length	20"   24.5"
Weight	.11 lbs   .16 lbs