

Economy Automatic Blood Pressure Monitor, Upper Arm

HCPCS: A4670-1,A4670



FIGURE A



FIGURE B

- WHO Classification
- Irregular Heart Beat Detection
- One Touch Operation
- Large Easy-to-Read Screen
- Date & Time Display
- Pulse Rate Indicator
- Digital Error Messages
- Travel Bag
- Auto-Off Function
- Memory Stores 120 Readings
- Average Last 3 Results

PRODUCTS

Description

BP3600	BP Monitor,Auto,Upper Arm,1/EA
--------	--------------------------------

SPECIFICATIONS

Description

Auto Power Off	3 Minutes After Last Reading
Carton Size	18.9" x 17.7" x 17.7"
Cuff Size	8.5"Ø14"
LCD Size	2.17" x 3.31"
Measurement Method	Oscillometric Method
Measurement Range	Pressure: 0mmHg~300 mmHg Pressure: +/- 3 mmHg
Memory Capacity	120 Readings with Date & Time
Operation Environment	50°F to 104°F, <85% RH
Power Supply	4x AAA Batteries or AC Adapter
Unit Dimensions	5.51" x 3.86" x 1.89"
Warranty	2 Years

ACCESSORIES

Description

BP2610	A/C Adapter, BP3400,BP3600,1/ea, EA
BP2000	Cuff,PlusSize,12.5"x17.75",BP3400,3600, EA
BP2001	Cuff,Std,BP3400,BP3600,BPHL,8.5-14, EA