

Foam Ring Cushion

HCPCS:



FIGURE A

- Provides excellent comfort and support while seated
- Open center is designed to provide maximum pressure relief
- Comes with durable, washable, zippered jacquard cover
- Made from soft, premium quality therapeutic foam
- Retail Packaged
- Dimensions: 16" (L) x 12 75" (W) x 2 75" (H)

PRODUCTS

Description

RTLPC23388	Invalid Ring - 16",1EA,6/CS
------------	-----------------------------

SPECIFICATIONS

Description

Material	Jacquard Knit Cover
----------	---------------------