

## Heating Pad



HCPCS:



FIGURE A



FIGURE B




FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Safe, Soothing & Fast-Acting Heat
- Absorbent sponge in pad provides moist heat that penetrates tissue quickly and safely
- Pad heats in seconds and monitors real-time temperature to prevent overheating
- Six temperature settings (113/122/131/140/149/158) and countdown timer settings (0:15/0:30/0:45 /1:00/1:30/2:00)
- Digital display features on/off switch, temperature and timer control, and indicator light
- Machine-washable pad
- Weighs 1.5 lbs. with large dimensions of 12 x 24 x 0.4
- Comes with portable tote bag for storage or convenient on-the-go use and elastic fastening strap
- 
- 
- 

### PRODUCTS

### Description

RTL19G007GR	Digital Heating Pad, 1EA, 6/CV
-------------	--------------------------------

### SPECIFICATIONS

### Description

Color	Gray
Warranty	Limited Lifetime